






Tested disinfectants

The manufacturers of the following cleaning and disinfecting agents have carried out material suitability tests and confirmed compatibility with our brush strips. In some cases, the tests were carried out in high doses to simulate compatibility over a long period of time. All disinfectants should be dosed and used according to the manufacturer's recommendations.

Manufacturer	Tested disinfectants	Tested concentration in %	Duration of the test	Approval by manufacturer	short-term use	long-term use
	BTS 2000	1,8	90 days			
	BTS 2001					
	Des AM blau					
	Virodex					
	Calgonit DS 683	5	19 days			
	Suma San con. D10.1 conc	2	30 days			
	SaniPerfect VS61	2	70 days			
	Suredis VT1	3				
	P3 topax 66	5	28 days			
	P3 topax 990	3				
	Deptal AS	2	33 days			
		20				
	1+1 Wofasteril SC super	1	30 days			
	MENNO® Clean	6	120 days			
	MENNO® Florades	8				
	A-QUASAN®					
	Lerades Clp 601	2	30 days			
	Lerasept AC	2				
	Lerasept L420	2,5				
	Lerasept SOC	6				

Tested disinfectants

The manufacturers of the following cleaning and disinfecting agents have carried out material suitability tests and confirmed compatibility with our brush strips. In some cases, the tests were carried out in high doses to simulate compatibility over a long period of time. All disinfectants should be dosed and used according to the manufacturer's recommendations.

Manufacturer	Tested disinfectants	Tested concentration in %	Duration of the test	Approval by manufacturer	short-term use	long-term use
	Blaue Sandy	100	4 days			
	Goldene Sandy					
	ASiRAL SRC	100	10 days	